




Isotope Batch Extraction System

IBEX Methane

Introduction


The IBEX methane is a fully automated system for the preparation and purification of CH₄ gas samples for clumped isotope analysis.

The IBEX can be integrated with most modern IRMS systems, including the Thermo Ultra. Gas connections are made directly to the IRMS bellows and/or systems inlet block.




How the IBEX Methane works


The IBEX uses a series of variable temperature silica traps and a GC-GC purification processes to clean CH₄ from gaseous samples.



The cryo-traps use a variable head-pressure system to raise and lower the liquid nitrogen level without physically moving the trap. This provides safe and accurate trap cooling and heating.



As the freezing point of CH₄ is close to LN₂ temperatures, a silica stationary phase is added to the IBEX cryo traps. This ensures near 100% recovery of CH₄ samples.





IBEX Control Centre Software

The IBEX Control Centre software (ICC) provides a fully integrated control and monitoring platform. Seamlessly linking with ISODAT sequence files, methods and scripting language.

Real-time control and readbacks of all functional objects are accessible in the main ICC screen, including valve state, trap temperatures and vacuum pressures.


The addition of custom controls and objects is also possible using the systems plug and play options. This includes pressure gauges, valves and other bespoke equipment.



Custom analytical methods can be created using the ICC method editor. This requires no prior scripting or coding knowledge. New methods can be created and deployed in minutes.

Existing methods can be edited with a few mouse clicks.

Advanced users can also edit and deploy methods using any XML editor.



IBEX functionality and specifications

CRYO TRAP PERFORMANCE
<p>The system includes three silica cryo-traps (two water traps and one cold finger).</p> <p>Water traps Max Temp +250 °C Min Temp -197 °C CO₂ release temperature -105 to -80°C Temperature increments 1°C</p> <p>IRMS cold finger Max Temp +130 °C Min Temp -197 °C CH₄ release temperature: +60°C Volume 200 ml (can be altered on request).</p>
LIQUID NITROGEN CONSUMPTION
<p>Approximately 45 litres per day.</p> <p>It is recommended that the customer provide a 90-litre (or larger) LN₂ storage vessel to supply the IBEX preparation system.</p>
SULPHUR TRAP
<p>Silver wool trap Volume 3 ml</p>
IBEX CONTROL CENTRE SOFTWARE
<p>Compatible with ISODAT Requires Windows 7 or above</p>
VACUUM SYSTEM
<p>1 x DX85 Edwards turbo molecular pump 2 x RV5 Edwards rotary pumps 1 x Edwards Penning gauge 2 x Edwards Pirani gauge 1 x Edwards TIC controller</p> <p>Please note alternative components may be used if they provide better performance than those described above.</p>
Vacuum Fittings
<p>All Vacuum joints and seals use VCR connections. All tubing is Ultron finished or equivalent</p>



Protium MS Ltd.
The Innovation Centre
Sci-Tech Daresbury
Keckwick Lane Daresbury
WA 4 4FS
www.protiumms.com



FS618518